

Percutaneous transcatheter valve replacement in individuals with mitral regurgitation unsuitable for surgery or transcatheter edge-to-edge repair: a prospective, multicountry, single-arm trial




Objectives

The ENCIRCLE trial is a prospective, single-arm, multi-center, pivotal trial to establish the safety and effectiveness of the SAPIEN M3 system in subjects with symptomatic MR $\geq 3+$ for whom commercially available surgical or transcatheter treatment options are deemed unsuitable due to clinical, anatomical, or technical considerations.¹

Methods

- 56 sites in the United States, Canada, Israel, the Netherlands, the United Kingdom, and Australia²
- Primary endpoint was a non-hierarchical composite of all-cause mortality or HF hospitalization at 1 year¹
- Secondary endpoints included the following at 1 year compared to baseline¹:
 - Improvement in MR
 - Improvement in NYHA functional class
 - Improvement in KCCQ-OS score
 - Decrease in LVEDVi



Key outcomes

- The **primary endpoint** of a non-hierarchical composite of all-cause mortality or heart failure hospitalization was achieved at 1 year with results **significantly better than the performance goal[†]** (25.2% vs 45%)²
- **<1% all-cause mortality** at 30 days²
- **96% of patients achieved MR $\leq 1+$** at 1 year²
- Patients treated with the SAPIEN M3 system had **significant improvements in health status** at 1 year²

¹Refer to Clinical Study Protocol for full enrollment criteria.

[†]Performance goal is based on the medical therapy arms of two trials with a similar subject population (MITRA-FR, COAPT).




Baseline characteristics

Baseline Characteristics ²	N=299
Age (years)	77.0
Male	51%
STS score, mitral valve replacement	6.6%
MV mean gradient (mmHg)	3.5
LVEF	49.5%
Congestive heart failure	75%
Prior CABG	30%
Prior TIA or stroke	19%
Prior mitral repair	9%
PPM or ICD	36%
Hypertension	84%
Atrial fibrillation	70%
NYHA class III/IV	71%
MR etiology	
FMR	58%
DMR	35%
Mixed	6%

Complex patient population at baseline with multiple comorbidities²

- 75% of patients had congestive heart failure and 70% had atrial fibrillation
- Mild to moderate MAC present in 24% of the Main Cohort patients
- Over 25% of patients had >1 reason for TEER unsuitability

Results



Procedural Information ²	SAPIEN M3 system N=299
Procedure time (min) ^{†,‡}	127.0 ± 47.1
Device time (min) [†]	102.8 ± 42.6
Dock deployment time (min) ^{**}	65.9 ± 35.0
Valve implant time (min) ^{††}	12.1 ± 14.3
Procedure aborted	4.0%
Conversion to surgery	0.0%
PVL closure	5.0%
ASD closure	
Routine closure	12.4%
Clinically significant closure	5.0%
Discharged home	96.3%
Index hospital stay	2.0 days

*Encircling difficulties, n=5; Low dock, n=5; Mobile echodensity, n=1; Unable to gain medial commissure access, n=1.

[†]Data on file. Sample sizes for procedure time measurements varied across groups.

[‡]Defined as the time from femoral vein access puncture to guide sheath removal.


[§]Defined as the time from guide sheath insertion to removal.

^{**}Defined as the time from dock delivery system insertion to removal.

^{††}Defined as the time from Commander delivery system insertion to removal.


Significant and sustained MR reduction²

Figure 1: MR Severity



Favourable left ventricular remodelling, with a significant reduction in LVEDVi^{2,3}


Figure 2: LVEDVi by Core Lab*



*By Core Lab, MedStar Health (Washington).

Achieved primary endpoint, a non-hierarchical composite of all-cause mortality or HFH, delivering results significantly better than the performance goal^{2*}

Figure 3: All-cause Mortality or HFH



Proven safety profile in comorbid, MR $\geq 3+$ patients²

Safety Outcomes [†]	30 Days	1 Year
All-cause mortality	0.7%	13.9%
• Cardiovascular mortality	0.7%	8.9%
Stroke	2.7%	9.3%
• Disabling stroke	1.7%	3.9%
• Non-disabling stroke	1.0%	5.5%
New Afib	7.9%	11.5%
New PPM	2.6%	5.5%
Major bleeding or above, MVARC [‡]	8.7%	18.5%
HF hospitalization	4.0%	16.7%
Mitral valve reintervention	2.3%	6.4%
Clinically significant valve thrombosis [§]	2.3%	6.7%
Valve embolization	0.0%	0.0%
Hemolysis requiring intervention	4.3%	7.1%

*Performance goal is based on the medical therapy arms of two trials with a similar subject population (MITRA-FR, COAPT).

[†]Data are Kaplan-Meier estimates %.


[‡]Major bleeding or above includes bleeding with MVARC primary bleeding scale of major, extensive, life-threatening or fatal.

[§]Clinically significant valve thrombosis includes leaflet thickening with impaired leaflet motion with mitral valve stenosis (increase in mean mitral valve gradient ≥ 5 mm Hg) and clinical signs or symptoms of mitral valve stenosis.

^{||}Hemolysis requiring blood transfusion or mitral reintervention.


Significant and sustained improvements in functional status²

Figure 4: NYHA Functional Class



Significant quality of life improvements²

Figure 5: KCCQ Overall Summary Score



Conclusion

The SAPIEN M3 system is an effective TMVR therapy with a proven safety profile for a wide range of patients²:




Significant and sustained MR reduction



Proven safety profile with low rates of all-cause mortality and HFH



Significant quality of life improvements



Abbreviations

Afib = Atrial fibrillation

ASD = Atrial septal defect

CABG = Coronary artery bypass graft

DMR = Degenerative mitral regurgitation

FMR = Functional mitral regurgitation

HF = Heart failure

HFH = Heart failure hospitalization

HR = Hazard ratio

ICD = Implantable cardioverter-defibrillator

KCCQ-OS = Kansas City Cardiomyopathy Questionnaire – Overall Summary

LVEDVi = Left ventricular end-diastolic volume index

LVEF = Left ventricular ejection fraction

MAC = Mitral annular calcification

MR = Mitral regurgitation

MV = Mitral valve

MVARC = Mitral Valve Academic Research Consortium

NYHA = New York Heart Association

PPM = Permanent pacemaker

PVL = Paravalvular leak

STS = Society of Thoracic Surgeons

TEER = Transcatheter edge-to-edge repair

TIA = Transient Ischemic Attack

TMVR = Transcatheter mitral valve replacement

1. The ENCIRCLE Trial. Clinicaltrials.gov Identifier: NCT04153292. Updated August 09, 2025. <https://clinicaltrials.gov/study/NCT04153292>
2. Guerrero M, et al. Percutaneous transcatheter valve replacement in individuals with mitral regurgitation unsuitable for surgery or transcatheter edge-to-edge repair: a prospective, multicountry, single-arm trial. *Lancet*. 2025 Nov 29;406(10519):2541-2550. doi: 10.1016/S0140-6736(25)02073-2
3. Daniels D. et al. Percutaneous Transcatheter Valve Replacement for Mitral Regurgitation: 1-Year Outcomes from the ENCIRCLE Trial. Presented at TCT 2025.

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, ENCIRCLE, the ENCIRCLE logo, SAPIEN, and SAPIEN M3 are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2026 Edwards Lifesciences Corporation. All rights reserved. PP--EU-11278 v2.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com



Edwards