




InteleViewer/ InteleConnect Installation and Setup



The following steps explain how to install and setup InteleViewer™ software. If you already installed InteleViewer™ , please proceed to Step 2.

Step 1: Install InteleViewer™

1. In your web browser, go to the following address:
<https://pacs.i-med.com.au/InteleViewer-tracks/enterprise/> (Note: address is CASE sensitive)
2. Download the relevant InteleViewer™ Installer (Windows 64-bit or Mac OS X)
For Windows: Choose Run from the popup window to install.
For Mac: Double-click the downloaded file to install and drop it into the applications folder.
3. Accept the Terms and Conditions and click Next.
4. Click Next on any further screens that pop up.
5. The Install Complete screen will appear when it has finished.



Step 2: Setup IntelViewer™

1. Click the IntelViewer™ icon on your desktop to open the Login Screen.
2. In the Login Screen, click the Edit Bookmarks button.
3. In the Server Bookmarks window, click the Add button.
4. In the New Bookmark window, enter the following server URL:

(VIC/NSW Metro areas, QLD, WA)

<https://pacs.i-med.com.au>

(VIC: Gippsland, Bass Coast, TAS, NT, NSW: Riverina Region, Albury/Wodonga, Wangaratta)

<https://pacs.ril.com.au>

(NSW inner west: Ashfield, Campsie, Drummoyne, Five Dock, Newtown, Newtown MRI, North Strathfield, Mater Imaging)

<https://www.iwpacs.com.au>

Enter Name: I-MED Radiology

- Do not hit ENTER / RETURN
 - Ensure the Protect Bookmark tick-box is checked
 - Check the Use Custom Server Name tick-box
5. Click Ok to close the New Bookmark window.
 6. Click Close to close the Server Bookmarks window.



Step 3: Login to IntelViewer™

1. In the Login Screen enter the username and password provided to you by I-MED Radiology. (To create a username please go to <https://portal.i-medonline.com.au/apply>).
2. If you are logging in for the first time you will receive an IntelBrowser prompt to accept the T&C's of the confidentiality agreement. Once accepted, close and reopen IntelViewer, then Login.

You are now ready to use the I-MED Radiology IntelViewer™ software.


If you would like to access images and reports without using the software, please see the following page for further information.



IntelConnect Online Access

Access your images and reports through a web browser on any computer without installing the IntelViewer™ software.

1. In your web browser, go to the following address:
<https://pacs.i-med.com.au/Portal/app#/>
2. Login with your username and password.
3. Search for your patient by 'surname/first name' and click the magnifying glass.




iPhone or iPad App*

Access your images and reports on your mobile devices through the IntelConnect™ app on your iPhone or iPad.

1. Download the IntelConnect™ app from the Apple App Store on your iPhone or iPad.
2. Open the app and enter the following details into the New Account screen:
Name: I-MED Radiology
Server: <https://pacs.i-med.com.au>
Username: your username
Password: your password
3. Click Done to add the account.
4. In the Login Screen enter your username and password.
5. Search for your patient by 'surname/first name' and click the magnifying glass.

***The IntelConnect™ app is currently only available for Apple devices.**



Support

Contact the I-MED Radiology Referrer Success Team if you require help setting up IntelViewer™, have not been provided with a username or password or need assistance with how to use IntelViewer™.

Phone: 1300 147 852 (during business hours)

Email: referrer.success@i-med.com.au