

ZenPure

PureFlo Nylon Screen Filter Capsules (Nylon Screen on PP Construction)

Screen Filtration

PureFlo® nylon screen capsule filters are designed for filtration and clarification applications in the pharmaceutical, biotechnology, food and beverage, medical, chemical, and water worldwide. This family of products is particularly suitable for the fluids with large sizes particles. The nylon screen removes particles by a sieving. Hard particles can be easily removed by the use of these units.

The nylon screen in an all-polypropylene construction provides excellent chemical compatibility. It can also be constructed in Nylon, Polyethylene and Gamma stabilized PP shells. PureFlo nylon screen filter capsules are well-suited for applications. No adhesives, binders, or surfactants are used in the manufacturing process. All capsules are rinsed with high-purity water to reduce extractables and downtime in a clean room environment to ensure filter performance each and every time out of the package.



Applications							
Clarification	Water						
Hard Particle	Beverages						
Cell Removal	Pharmaceuticals						
Chemicals	Biologics						

Specifications

Materials of Construction:	Media: Nylon Screen Media Supports: None Cage, Core, End Caps: Polypropylene, Gamma Stabilized PP, Nylon Sealing: Thermally bonded
Fitting Connections:	See ordering guide for the availability. (Custom adaptors available upon request)
Effective Filtration Area:	Lengths: 1.5" 0.65ft² (500cm²), 2.5" 1.3ft² (1000cm²), 5" 2.1ft² (2000cm²), 7.5" 3.2ft² (3000cm²) and 10" 4.3ft² (4000cm²) Diameter: 2.75 in. (70mm)
Available Ratings:	10um, 20um, 40um, 60um, 100um, and 200um
Operating Conditions:	Maximum Operating Pressure: Liquid: 5.5 bar (80psi) at 77°F/25°C Gas: 4.1 bar (60psi) at 77°F/25°C Maximum Forward Differential Pressure: 4.21 bar (60psi) at 68°F/20°C Maximum Reverse Differential Pressure: 2.1 bar (30 psi) at 68°F/20°C Minimum Burst Pressure: 8.3 bar (120psi) at 77°F/25°C Maximum Operating Temperature: 176°F (80°C)
Regulatory Compliance:	The filters are constructed with polypropylene resins and filtration media in compliance with 21CFR Part 177 of the US Code of Federal Regulations and USP Class VI Biological Test for Plastic.


ZenPure

PureFlo Nylon Screen Filter Capsules

Specifications (cont.)

Sanitization:

The filters can be sterilized by autoclaving for up to 25 cycles at 125°C (250°F) for 60 minutes. The filters can also sanitized by hot water or common chemicals that are compatible with Nylon and polypropylene.



60um Nylon Screen

PureFlo Nylon Screen Filter Capsule Ordering Guide

PureFlo Capsule Filters	Vent & Drain Fitting	Filter Media	Pore size (Micron)	Length	Input Fitting	Output Fitting	Options		
	L = Female Luer		100 = 10	H = 1.5"	1H = 1/8" Hose barbs	1H = 1/8" Hose barbs	-GP = Gamma stable	O-Ring for Quick	
70mm	Lock		200 = 20	S = 2.5"	1NF = 1/8" FNPT	1NF = 1/8" FNPT	polypropylene shell	Coupling	
	P = M8 Threaded	NS = Nylon Screen	400 = 40	L = 5"	1Q = 1/8" Male Quick Coupling for Metal latch	1Q = 1/8" Male Quick Coupling for Metal latch	-NY = Nylon Construction	-OE = O-ring EPDM	
	Plug - <u>Limited</u>	Screen	600 = 60	E= 7.5"	1QF = 1/8" Female Quick Coupling w/Metal latch	1QF = 1/8" Female Quick Coupling w/Metal latch	-FB = Filling bell	-ON = O-ring Nitrile	
	<u>options</u>		10X = 100	K = 10"	2C = 1/4" Compression	2C = 1/4" Compression	-FC = Filling bell with Cap-	DV = O-ring Viton	
Capsule	N = No vent or		20X = 200		2H = 1/4" to 1/2" Hose barbs	2H = 1/4" to 1/2" Hose barbs	-ETO = Ethylene oxide	Blank = O-ring Silicon	
	drain fittings				2HS = 1/4" to 3/8" Hose barbs	2HS = 1/4" to 3/8" Hose barbs	sterilization	(Standard)	
					2LF = Large Female Luer Lock fitting	2LF = Large Female Luer Lock fitting	-NI = No vent or drain	PreFilters	
					2N = 1/4" MNPT	2N = 1/4" MNPT	inlet fitting	Prefliters	
					2NF = 1/4" FNPT	2NF = 1/4" FNPT	-NO = NO vent or drain	-G(pore Size) = Glass	
					2Q = 1/4" Male Quick Coupling for Metal latch	2Q = 1/4" Male Quick Coupling for Metal latch	Outlet fitting	Fiber PreFilter	
					2QF = 1/4" Female Quick Coupling w/Metal latch	2QF = 1/4" Female Quick Coupling w/Metal latch	-E = Polyethylene shell and HDPE media support for gamma sterilization	-P(pore Size) = PolyPro Media PreFilter	
					2QFP = 1/4" Female Quick Coupling w/Plastic latch	2QFP = 1/4" Female Quick Coupling w/Plastic latch			
					2QP = 1/4" Male Quick Coupling for Plastic latch	2QP = 1/4" Male Quick Coupling for Plastic latch		-S(pore Size) = PES	
					3C = 3/8" Compression	3C = 3/8" Compression		PreFilter	
					3H = 3/8" Hose Barb	3H = 3/8" Hose Barb			
					3HE = 3/8" Hose Barb Elbow	3HE = 3/8" Hose Barb Elbow			
					3NF = 3/8" FNPT	3NF = 3/8" FNPT			
					4C = 1/2" Compression	4C = 1/2" Compression			
					4H = 1/2" Hose Barb	4H = 1/2" Hose Barb			
					4N = 1/2" MNPT	4N = 1/2" MNPT			
					4QP = 1/2" Male Quick Coupling for Plastic latch	4QP = 1/2" Male Quick Coupling for Plastic latch			
					MT = 1/2" Tri clamps	MT = 1/2" Tri clamps			
					TC = 1-1/2" Tri clamp	TC = 1-1/2" Tri clamp			
Example - SK Series Capsule, Luer Lock V/D, 2.5" Internal filter, 20 micron Nylon Screen Capsule, Standard Grade, with 1.5" Tri-Clamp fittings would be SKLNS200STCTC									
ick Couplings are molded. Female Quick Coupling with Plastic latch will be threaded on, compatible with CPC (Colder), LinkTech and others.									

uick Couplings are molded, Female Quick Coupling with Plastic latch will be threaded on, compatible with CPC (Colder), LinkTech and others.



www.internationalfilterproducts.com

7914 Ajay Dr. Sun Valley, CA 91352 Tel: (818) 504-8115 - Office Tel: (818) 917-9538 - Cell

E-mail: sales@internationalfilterproducts.com

 $ZenPure\ and\ PureFlo\ are\ registered\ trademarks\ of\ ZenPure\ Corporation\ or\ an\ affiliated\ company.$ Copyright 2003-2012\ ZenPure\ or\ an\ affiliated\ company. All rights reserved.



