# DEPTH-CLEAR<sup>™</sup> P & LE Series

Depth Filter Modules for Biopharmaceutical Applications



#### Description

DEPTH-CLEAR P Series depth filter modules are specifically formulated to meet the performance demands of biopharmaceutical and biotech applications. The filter media utilize mechanical and electrokinetic adsorptive capture mechanisms to remove particles, microorganisms, colloids, and pyrogen from critical process streams. Purolator, a leader in depth media development, has invested years of research and development to create a range of filter grades to meet stringent application requirements.

The DEPTH-CLEAR proprietary formulation process produces filter media with exceptionally high, effective filtration area that provides superior particle retention, high contaminant holding capacity, and long on-stream filter life. Manufacturing by means of an advanced, highly automated production process results in very consistent product quality and filtration performance. Purolator can provide technical assistance to validate the performance of DEPTH-CLEAR filters in new processes and applications.

DEPTH-CLEAR Pharmaceutical Grade (P Series) filter modules are manufactured to procedures described in a Drug Master File (DMF) on record at the National Center for Drugs and Biologics. The filter media is non-toxic in accordance with systemic injection testing. A special low extractable (LE Series) grade is available for susceptible processes and high purity applications. The DEPTH-CLEAR P & LE grades have complete product control and traceability, as set forth in the Drug Master File.

#### Benefits


- High surface area filter modules provide superior filtration performance and long on-stream life cycles.
- Proprietary formulation process results in depth media with exceptional particle retention via mechanical and electrokinetic capture.
- The advanced manufacturing process control results in very uniform filtration media and consistent lot-to-lot performance.
- Non toxic Meets the requirements of USP XXIII Class VI Biological Tests for Plastics.
- Special low extractable formulation (LE Series) available for sensitive applications.
- Special formulations available for high temperature applications (>140°F).
- Filter modules available in a wide variety of dimensions and configurations to fit most standard housings.

### Applications

- Antibiotics
- Blood Products
- Intravenous Solutions
- Cell Separation
- Cough Syrups
- Chromatography Column Protection
- Post Fermentation Clarification
- Small and Large Volume Parentals



#### Performance Specifications



Differential Pressure to 35 psid (2.5 bar) @ 140° F (60° C)

Maximum Operating<br/>TemperatureR Grade: 140° F (60° C)K Grade: 185° F (85° C)Recommended0.5 to 1 GPM/1ft²Flow Rate20 to 40 LPM/1m²

FDA Conformity

All materials conform to FDA standards regarding material contact during food

and beverage processing.

Toxicity DEPTH-CLEAR Series filter components

comply with USP XXIII standards for

Class VI-121 plastics.

#### Materials of Construction

Media Cellulose Fibers with Inorganic Filter Aids

End Caps Polypropylene
Cartridge Straps Steel


Gaskets/O-Rings Ethylene Propylene, Silicone, Neoprene, Nitrile, Viton®, Teflon®


Dimensions (nominal)

| Grade Dimensions         | NAXXXXC9-XX     | NAXXXXC16-XX   | NAXXXXD14-XX   |  |
|--------------------------|-----------------|----------------|----------------|--|
| DIAMETER<br>inches<br>cm | 11 1/8<br>28.3  | 11 1/8<br>28.3 | 16 3/4<br>42.5 |  |
| HEIGHT inches cm         | 7 11/16<br>19.5 | 10 7/8<br>27.6 | 10 7/8<br>27.6 |  |
| FILTER AREA<br>ft²<br>m² | 11<br>1.0       | 19<br>1.8      | 38<br>3.5      |  |
| # CELLS                  | 9               | 16             | 14             |  |

Viton and Teflon are registered trademarks of E.I. duPont de Nemours and Company. Depth-Clear and The Clear Advantage are trademarks of Clarcor, Inc.

#### Flow Rate vs Differential Pressure





The above graphs indicate pressure drop versus flow rate of the various grades of Depth-Clear filter media in water. Adjustments may be required for other fluids. Please check with your Purolator representative for assistance in sizing.

**Differential Pressure** 

## Cartridge Ordering Information

| Media Type   | Grade Designation                                 | Formulation        | Series                                      | Handles                               | # of Cells         | Gasket Mate   | erial                                  |
|--------------|---|--------------------|---|---------------------------------------|--------------------|---|--|
| NA           | 10  | K                  | Р   | 05                                    | C9                 | -01   |  |
| Non-Asbestos | 10, 15, 30, 45, 60, 90, 150, 300, 600, 1000, 2000 | R Grade<br>K Grade | P = Pharmaceutical<br>LE = Low Extractables | 05 = Handles<br>No Symbol = No Handle | C9<br>esC16<br>D14 | EPR<br>Neoprene<br>Silicone<br>Nitrile<br>Viton<br>Teflon | -01<br>-02<br>-03<br>-04<br>-08<br>-09 |



7914 Ajay Dr. Sun Valley Dr., CA 91352 www.international lterproducts.com sales@international lterproducts.com Phone: (818) 504-8115

Fax: (818) 841-2947