

HemoSphere Advanced Monitor ClearSight Module

Competency checklist

Name: _____

Unit: _____

Date: _____


Please visit:



Edwards
Lifesciences
Clinical
Education



YouTube
channel
for video
instructions



Critical elements

ClearSight module setup

	Met	Not met
Plug the HemoSphere power cord to a hospital grade power outlet.	<input type="checkbox"/>	<input type="checkbox"/>
Ensure the HemoSphere ClearSight module is present and connected to the HemoSphere advanced monitor.	<input type="checkbox"/>	<input type="checkbox"/>
Press the power button on the front of the monitor to power on the HemoSphere advanced monitor.	<input type="checkbox"/>	<input type="checkbox"/>
Once prompted, select either the "New Patient" button and enter patient data or the "Continue Same Patient" button, then select "Next."	<input type="checkbox"/>	<input type="checkbox"/>
Plug the pressure controller into the HemoSphere ClearSight module.	<input type="checkbox"/>	<input type="checkbox"/>
Wrap the pressure controller band around the patient's wrist (either wrist can be used however the non-dominant arm is preferred).	<input type="checkbox"/>	<input type="checkbox"/>
Attach the pressure controller to the band.	<input type="checkbox"/>	<input type="checkbox"/>
Select the proper size finger cuff by using the ClearSight finger sizing aid, best results are obtained from the middle, ring, or index finger.	<input type="checkbox"/>	<input type="checkbox"/>
Open the finger cuff and place the finger cuff between the lines, with the finger cuff centered between the knuckles.	<input type="checkbox"/>	<input type="checkbox"/>



Edwards

ClearSight module setup (continued)	Met	Not met
Gently lead the cuff between the two fingers.	<input type="checkbox"/>	<input type="checkbox"/>
Tightly wrap the cuff around the finger; do not rotate the cuff after application.	<input type="checkbox"/>	<input type="checkbox"/>
Connect the finger cuff to the pressure controller.	<input type="checkbox"/>	<input type="checkbox"/>
Connect the heart reference sensor (HRS) to the pressure controller.	<input type="checkbox"/>	<input type="checkbox"/>
Vertically align both ends of the HRS and touch the on-screen zero button. Wait for indication that the HRS has been zeroed.	<input type="checkbox"/>	<input type="checkbox"/>
Attach the end of the HRS with a heart to the patient at the phlebostatic axis level by using an HRS clip.	<input type="checkbox"/>	<input type="checkbox"/>
Attach the other end of the HRS to the finger cuff.	<input type="checkbox"/>	<input type="checkbox"/>
Touch the start monitoring icon on the navigation bar to begin monitoring.	<input type="checkbox"/>	<input type="checkbox"/>
Touch the home icon the return to the main monitoring screen.	<input type="checkbox"/>	<input type="checkbox"/>
Check the status bar to ensure there are no faults or alters that may indicate a problem in establishing the arterial waveform. If a displayed message appears, touch the question mark icon for more information on the displayed message.	<input type="checkbox"/>	<input type="checkbox"/>
To select desired monitoring screen view, touch the settings icon and choose "Select Screens" tab.	<input type="checkbox"/>	<input type="checkbox"/>
To change the parameter tiles shown on the right side of the screen, touch the screen area where the parameter label is displayed. The tile configuration menu will show the selected parameter highlighted in color. Available parameters appear on the screen without highlights; touch any parameter to select the replacement parameter.	<input type="checkbox"/>	<input type="checkbox"/>
To adjust parameter alarms/targets, touch the settings icon, then select advanced setup (Contact your hospital administrator or IT department for passwords), then select parameter settings, and finally select alarms/targets. Once completed, select the home icon to apply settings and return to the main monitoring screen.	<input type="checkbox"/>	<input type="checkbox"/>

Passed competency

Needs to repeat competency

Validated by: _____

Date: _____

CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards, Edwards Lifesciences, the stylized E logo, ClearSight, and HemoSphere are trademarks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2021 Edwards Lifesciences Corporation. All rights reserved. PP--US-5658 v1.0

Edwards Lifesciences • One Edwards Way, Irvine CA 92614 USA • edwards.com



Edwards