Clinical summary

Transcatheter valve replacement with the EVOQUE system in patients with symptomatic ≥ severe tricuspid regurgitation




Objectives

The TRISCEND II trial is a prospective, multi-center, randomized pivotal trial evaluating the safety and effectiveness of transcatheter tricuspid valve replacement using the Edwards EVOQUE system with optimal medical therapy compared to optimal medical therapy alone in patients with symptomatic ≥ severe tricupsid regurgitation (TR).

Methods



- TR graded as ≥ severe
- Symptomatic despite optimal medical therapy



- 45 enrolling sites in the United States and Germany
- Primary safety and effectiveness endpoint was a hierarchical composite assessed at 1 year:
 - All-cause mortality,
 - RVAD implantation or heart transplant,
 - Tricuspid valve surgery or percutaneous tricuspid intervention,
 - Annualized rate of HF hospitalizations,
 - − KCCQ-OS improvement $\Delta \ge 10$ points,
 - − NYHA improvement $\Delta \ge 1$ functional class, and
 - − 6MWD improvement $\Delta \ge 30$ meters

* 12.6±2.1% for EVOQUE TTVR vs 15.2±3.3% for medical therapy alone.
 † 20.9±2.6% for EVOQUE TTVR vs 26.1±4.1% for medical therapy alone.

Key Outcomes

- The primary endpoint was met, demonstrating
 EVOQUE TTVR was superior to medical therapy alone
 with a win ratio of 2.02
 (95% Cl, 1.56 to 2.62; p<0.001).
- 95.3% achieved TR reduction to ≤mild with EVOQUE TTVR compared to 2.3% with medical therapy alone at 1 year.
- There were favorable numerical trends for
 EVOQUE TTVR compared to medical therapy alone in all-cause mortality (12.6% vs. 15.2%)* and heart failure hospitalizations (20.9% vs 26.1%)[†] at 1 year.
- Patients treated with
 EVOQUE TTVR had greater
 improvements in health
 status, including 18.4-point
 increase in KCCQ-OS,
 91.0% in NYHA class I/II, and
 23.2m increase in 6MWD
 at 1 year.



Results

Baseline Characteristics	EVOQUE TTVR n=259*	Medical therapy alone n=133
Age, years	79.3	79.1
Female	74.9%	76.7%
NYHA class III-IV	73.0%	69.2%
Atrial fibrillation	96.1%	92.5%
Chronic kidney disease	54.1%	59.4%
Pacemaker/ICD	38.2%	39.8%
Ascites	18.5%	21.8%
KCCQ-OS (mean)	52.8	50.6
Secondary TR etiology	74.1%	71.4%

Procedural Outcomes

- 100% performed via percutaneous femoral vein access
- No procedural deaths and 1.2% conversion to surgery
- 93% discharged home


Procedural Characteristics	n=259*
Device implanted as intended	95.4% [†]
Device time, median (IQR), min	56.5 (41.0, 75.0)
Length of stay, median (IQR), days	3.0 (2.0, 6.0)


* Of the 267 patients randomized to EVOQUE TTVR, 8 exited prior to procedure

⁺ 12 patients (4.6%) did not receive device due to challenging imaging or anatomy


Echocardiographic Outcomes

■ EVOQUE TTVR achieved 95.3% ≤mild TR compared to 2.3% in medical therapy alone. TR was eliminated in 72.6% of patients who received EVOQUE TTVR.






Results (continued)



Results (continued)

Primary Endpoint

Based on 34447 possible pairs, there were 21397 wins for EVOQUE TTVR, 10591 wins for medical therapy alone, and 2459 ties, resulting in a win ratio of 2.02 (95% CI, 1.56 to 2.62; p<0.001). The primary safety and effectiveness endpoint was met, demonstrating EVOQUE TTVR was superior to medical therapy alone.</p>





Results (continued)

30-day Safety Outcomes

CEC-Adjudicated Safety Events	EVOQUE TTVR n=259
Cardiovascular mortality	3.1%
Myocardial infarction	0.8%
Stroke	0.4%
Severe bleeding*	10.4%
Non-elective tricuspid re-intervention	0.8%
New pacemaker/CIED implants in pacemaker-naïve patients	24.7%

* Fatal, life-threatening, extensive, or major bleeding, as defined by Mitral Valve Academic Research Consortium (MVARC; Stone et al. 2015).

Conclusion

- Patients with severe TR experience significant symptom burden with diminished quality of life.
- TRISCEND II is the first randomized controlled trial studying tricuspid valve replacement compared to medical therapy alone.
- Results from the TRISCEND II trial establish TTVR as an effective therapy with a proven safety profile for patients with symptomatic
 ≥ severe TR, with consistent TR resolution accompanied by meaningful health status benefits.

Abbreviations (alphabetical)

- 6MWD = Six-minute walk distance CEC = Clinical events committee CI = Confidence interval CIED = Cardiac implantable electronic device HF = Heart failure ICD = Implantable cardioverter defibrillator IQR = Interquartile range IVC = Inferior vena cava KCCQ-OS = Kansas City Cardiomyopathy Questionnaire – Overall Summary LVOT = Left ventricular outflow tract MVARC = Mitral Valve Academic Research Consortium NYHA = New York Heart Association RVAD = Right ventricular assist device RVEDD = Right ventricular end-diastolic diameter
- **SV** = Stroke volume
- **TR** = Tricuspid regurgitation
- **TTVR** = Transcatheter tricuspid valve replacement





Experience REVOLUTIONARY

Learn more at Edwards.com/gb/EVOQUE

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, Edwards EVOQUE, EVOQUE, TRISCEND II, and the TRISCEND II logo are trademarks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

 ${\rm @}$ 2024 Edwards Lifesciences Corporation. All rights reserved. PP--EU-9275 v1.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com

