

One-year outcomes of the TriCLASP post-market study


At one-year the PASCAL system demonstrated significant and sustained clinical benefits for patients with clinically significant TR:

- ✓ High Survival and low rates of HF hospitalisation
- ✓ Significant and sustained TR reduction
- ✓ Clinically meaningful improvements in functional status and quality-of-life
- ✓ Low rates of MAEs

Study design

Prospective, multicentre, single-arm, post-market clinical follow-up study to evaluate the **safety and effectiveness of the PASCAL system** in patients with **symptomatic, severe or greater TR**, in a European post-market setting.

Patient enrollment



Baseline characteristics

300 patients, 80±6 years old	78% NYHA III/IV
81% FTR / mixed ^b	91% Atrial fibrillation
76% TR≥ severe ^c	62% Kidney disease

^a1 due to complex anatomy, 1 had cardiac arrest during the procedure and the device was not implanted because grasping was deemed not possible.

^bPrimary 15.7%, pacer related 2.0%, indeterminate 1.3%. ^cCore laboratory: Cardialysis, Rotterdam, The Netherlands.


Low Major Adverse Events rate at 30 days¹

1.7% Composite MAE rate	1.3% All-cause mortality	2.3% SLDA rate
-----------------------------------	------------------------------------	--------------------------

¹n=291 patients, denominator includes patients who had an event and/or were followed to at least 30 days. Patients may have had more than one event.

Significant and sustained TR reduction to 1 year


87%
of patients
with TR≤2+
at 1-year with the
PASCAL system



Graph shows unpaired data.


Low rate of all-cause mortality and HF hospitalisation with significant reduction in annualised rate of HF hospitalisation to 1 year

Freedom from all-cause mortality and HF hospitalisation



88% Freedom from all cause mortality at 1-year


Annualised rate of HF hospitalisation



72% Relative reduction in annualised HF hospitalisation rate


Significant improvements in functional and quality of life outcomes at 1 year

NYHA Class



75% of patients had achieved NYHA Class I/II at 1-year with the PASCAL system

KCCQ Score



+8.3 Point improvement on KCCQ from baseline to 1 year

Conclusion



One-year outcomes from the TriCLASP study confirm sustained safety and effectiveness of the PASCAL system in patients with clinically significant TR in a post-market setting.

TR: tricuspid regurgitation; FTR: functional tricuspid regurgitation; MAE: major adverse events; HF: Heart failure; CEC: clinical events committee; SLDA: single leaflet device attachment; KCCQ: Kansas City Cardiomyopathy Questionnaire; NYHA: New York Heart Association

For details on statistical analyses please see reference 1

Reference

1. J. Hausleiter. Transcatheter Tricuspid Valve Repair: TriCLASP Study 1-Year Results. PCRLV 2024, Tricuspid Hotline. 25th Nov. 2024.

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, and PASCAL are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2025 Edwards Lifesciences Corporation. All rights reserved. PP--EU-9645 v1.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com



Edwards