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
Références

1. Davies, et al. Ability of an arterial waveform analysis-derived hypotension prediction index to predict future hypotensive events in surgical patients. *Anesthesia & Analgesia* 2019; 130(2) : 352-359.
2. Wijnberge, et al. Effect of a machine learning-derived early warning system for intraoperative hypotension vs standard care on depth and duration of intraoperative hypotension during elective noncardiac surgery: *The HYPE Randomized Clinical Trial*. *JAMA* 2020.
3. Schneck, et al. Hypotension prediction index based protocolized haemodynamic management reduces the incidence and duration of intraoperative hypotension in primary total hip arthroplasty: a single center feasibility randomized blinded prospective interventional trial. *J Clin Monit Comput* 2020.
4. Hatib, et al. Machine-learning algorithm to predict hypotension based on high-fidelity arterial pressure waveform analysis. *Anesthesiology* 2018; 129 : 663-674.
5. Ostadal, et al. Continual measurement of arterial dP/dt_{max} enables minimally invasive monitoring of left ventricular contractility in patients with acute heart failure. *Crit Care* 2019; 23:364
6. De Hert, et al. Evaluation of Left Ventricular Function in Anesthetized Patients Using Femoral Artery dP/dt_{max} . *Journal of Cardiothoracic and Vascular Anesthesia* 2006; 20(3): 325-330.
7. Tartiere, et al. Non-invasive radial pulse wave assessment for the evaluation of left ventricular systolic performance in heart failure. *Eur Journal of Heart Failure* 2007; 9: 477-483.
8. Garcia Monge, et al. Performance comparison of ventricular and arterial dP/dt_{max} for assessing left ventricular systolic function during different experimental loading and contractile conditions. *Crit Care* 2018.
9. Garcia Monge, et al. Dynamic arterial elastance as a predictor of arterial pressure response to fluid administration: a validation study. *Crit Care* 2014; 18: 626-637.
10. Cecconi, et al. The use of Pulse Pressure Variation and Stroke Volume Variation in spontaneously breathing patients in spontaneously breathing patients to assess dynamic arterial elastance and to predict arterial pressure response to fluid administration. *Anest & Analg* 2015; 120: 76-84.
11. Garcia Monge, et al. Dynamic arterial elastance to predict arterial pressure response to volume loading in preload-dependent patients. *Crit Care* 2011; 15: R15

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© 2020 Edwards Lifesciences Corporation. Tous droits réservés. PP--US-4928 v1.0