




## Objectives


The TRISCEND II pivotal trial met its primary endpoint, and results have previously been reported. The objective of the TRISCEND II trial is to compare the health status outcomes of patients with symptomatic  $\geq$  severe tricuspid regurgitation (TR) treated with transcatheter tricuspid valve replacement (TTVR) with the Edwards EVOQUE system plus optimal medical therapy compared to those who received optimal medical therapy alone.

<b>TRISCEND II Trial Design</b>		N=400 randomized (2:1) EVOQUE TTVR vs medical therapy alone	
<b>Health Status Assessment</b>		<b>Kansas City Cardiomyopathy Questionnaire (KCCQ)</b> at baseline, 30 days, 6 months, 1 year	
<b>Baseline Characteristics*</b>		<b>Mean Age</b> 79.2 years	<b>Sex</b> 75.5% female
		<b>KCCQ-OS</b> 52.1 points	<b>TR Grade</b> 56.1% massive/torrential

\*mITT (modified intent to treat safety population; n=392)

## Significantly improved KCCQ-OS at 30 days, with further improvements at 6 months that were sustained to 1 year

Figure 1. KCCQ-OS over 1 year




## Key Outcomes

- EVOQUE TTVR patients compared to those who received medical therapy alone experienced **significantly greater improvement in KCCQ-OS at 30 days** ( $\Delta 11.8$  points), with **further improvement through 6 months** ( $\Delta 20.8$ ) that was **sustained to 1 year** ( $\Delta 17.8$ ).
- Twice as many patients were "alive and well"** (KCCQ-OS  $\geq 60$  and no decline  $\geq 10$  points from baseline) and **"alive with large improvement"** ( $\Delta \geq 20$  points from baseline) at 1 year after EVOQUE TTVR than with medical therapy alone, which was statistically significant.
- Patients with  $\geq$ severe TR at baseline experienced significant improvement in KCCQ-OS at 1 year with EVOQUE TTVR compared to medical therapy alone, with the **greatest improvements seen in patients starting with massive and torrential TR** ( $\Delta 22.6$  and  $\Delta 23.3$ , respectively).
- Every 10-point increase in KCCQ-OS from baseline to 30 days was associated with a significant 19% reduction in the composite death or HF hospitalization risk at 1 year**, regardless of treatment group.




**Twice as many patients were “alive and well” and “alive with large improvement” at 1 year after EVOQUE TTVR than with medical therapy alone**

**Figure 2. Survival and health status at 1 year**




**A greater proportion of patients experienced large or moderate changes in KCCQ-OS with EVOQUE TTVR; conversely, a greater proportion of patients experienced worse outcomes on medical therapy alone**

**Figure 3. Change in health status by KCCQ-OS at 1 year**



**Patients with more severe TR, better exercise capacity (6MWD), and preserved right ventricular function (RV FAC) at baseline had significantly greater KCCQ-OS benefit at 1 year after EVOQUE TTVR**

**Figure 4.** KCCQ-OS sub-group analyses by baseline characteristics



**Every 10-point increase in KCCQ-OS from baseline to 30 days was associated with a significant 19% reduction in the composite death or heart failure hospitalization risk at 1 year, regardless of treatment group**


**Table 1.** Association of KCCQ-OS at 30 days with subsequent risk of death or HF hospitalization through 1 year

	Adjusted Hazard Ratio (95% CI)	p-value	Treatment interaction p-value
Death	0.87 (0.74-1.03)	0.100	0.327
Heart failure hospitalization	0.80 (0.71-0.90)	<0.001	0.419
Death or heart failure hospitalization	0.81 (0.73-0.90)	<0.001	0.945

Hazard ratios are scaled per 10-points on the KCCQ-OS. Adjusted models included age, sex, and chronic obstructive lung disease.

**Conclusion**

- Patients with symptomatic ≥ severe TR experience substantial impairment in health status.
- Compared with medical therapy alone, treatment of patients with symptomatic ≥ severe TR with EVOQUE TTVR resulted in significant and sustained improvements in patients’ symptoms, function, and quality of life.
- Significant health status benefits were evident at 30 days after EVOQUE TTVR, continued to increase through 6 months, and remained durable through 1 year.



**Experience REVOLUTIONARY**  
Learn more at [Edwards.com/gb/EVOQUE](https://Edwards.com/gb/EVOQUE)

**Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult [eifu.edwards.com](https://eifu.edwards.com) where applicable).**

Edwards, Edwards Lifesciences, the stylized E logo, Edwards EVOQUE, EVOQUE, TRISCEND II, and the TRISCEND II logo are trademarks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2024 Edwards Lifesciences Corporation. All rights reserved. PP--EU-9276 v1.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • [edwards.com](https://edwards.com)



**Edwards**