A photograph of a woman with short white hair and a young girl with blonde hair, both smiling and looking upwards. The woman is holding the girl. They are outdoors, with a blurred background of greenery and a warm sunset sky. A red rectangular overlay is positioned on the left side of the image, containing white text.

PASCAL System for the treatment of tricuspid regurgitation

Body of evidence

Clinical evidence for the use of the PASCAL system
for the treatment of tricuspid regurgitation



Edwards

PASCAL System body of evidence: Tricuspid regurgitation

Consistent favourable safety and effectiveness outcomes with the PASCAL system¹⁻⁴



2,000+ patients
targeted in studies



3 studies and **1** registry*



4+ years of study[†]

TR severity evaluated
by independent ECL
Prespecified adverse
events adjudicated
by CEC

THE CLASP TR STUDY



Prospective, multicentre, single-arm
early feasibility study¹

65
patients

1
year[†]

Significant and sustained improvements in TR, functional status and quality of life at 1 year.¹



Figure shows unpaired data. Adapted from Kodali KS et al. 2023. For details on statistical analyses, please see reference 1.

88%
freedom from
all-cause mortality

79%
freedom
from HFH

92%
NYHA class I/II

+18 pts
KCCQ-OS

THE CLASP II TR TRIAL



Prospective, multicentre, randomised,
controlled pivotal trial (roll-in cohort)²

73
patients[†]

30
days[†]

Significant TR reduction with no adverse safety signals, mortality or hospital readmissions in the roll-in cohort. The randomised trial is ongoing.²




Figure shows paired analysis. Adapted from Young M. 2022. For details on statistical analyses, please see reference 2.

100%
freedom from
all-cause mortality

100%
freedom
from HFH

86%
NYHA class I/II

+18 pts
KCCQ-OS

PASTE registry




Investigator-initiated, multicentre, retrospective and prospective, observational cohort study³

1,059
patients

1
year

Significant and sustained TR reduction and clinical improvements at 1 year in a real-world, high-risk patient population with complex anatomy at baseline:³

Mean TRI-SCORE risk **23%** | Coaptation gap >8 mm **24%** | Transvalvular CIED lead **27%** | >3 leaflets **43%**



86%
freedom from
all-cause mortality



84%
freedom
from HFH



66%
NYHA class I/II



-9 pts
MLHFQ

Figure shows unpaired analysis. Adapted from Wild MG et al. 2025. For details on statistical analyses, please see reference 3.

THE TriCLASP STUDY




Prospective, multicentre, single-arm, post-market clinical follow-up⁴

300
patients

1
year[†]

Significant and sustained improvements in TR reduction, functional status and quality of life at 1 year in patients with clinically significant TR in a post-market setting.⁴



88%
freedom from
all-cause mortality



83%
freedom
from HFH



75%
NYHA class I/II



+8 pts
KCCQ-OS

Figure shows unpaired analysis. Adapted from Hausleiter J. 2024. For details on statistical analyses, please see reference 4.

*CLASP TR, TriCLASP, CLASP II TR and PASTE registry.

[†]Based on reported follow-up.

Edwards PASCAL Precision System

Accurate, intuitive control

Advanced catheter and handle design facilitates smooth navigation and implant deployment*

Versatile implant configuration


Adapt to specific procedural and anatomical needs

Atraumatic clasp and closure

Enhance leaflet capture with atraumatic reclasp capabilities


Predictable release

Deploy the implant with procedural confidence†



*Design data on file and marketing evaluation.

†Performance and design data on file.



Learn more about the PASCAL system
at Edwards.com/gb/PASCAL

For details on statistical analyses, please see references 1–4.

Abbreviations

CEC: clinical events committee
CIED: cardiac implantable electronic device
ECL: echocardiography core lab
HFH: heart failure hospitalisation

KCCQ-OS: Kansas City Cardiomyopathy Questionnaire overall summary score

MLHFQ: Minnesota Living with Heart Failure Questionnaire
NYHA: New York Heart Association
TR: tricuspid regurgitation

References

1. Kodali SK, Hahn RT, Davidson CJ et al. 1-Year outcomes of transcatheter tricuspid valve repair. *J Am Coll Cardiol*. 2023; **81**: 1766–76.
2. Young MN. Transcatheter tricuspid valve repair: The CLASP II TR trial roll-in cohort. PCR London Valves, 2022, London, UK.
3. Wild MG, Stolz L, Rosch S et al. Transcatheter valve repair for tricuspid regurgitation: 1-year results from a large European real-world registry. *J Am Coll Cardiol*. 2024; 2025; **85**: 220–31.
4. Hausleiter J. Transcatheter tricuspid valve repair: TriCLASP study 1-year results. PCR London Valves, 24–26 November 2024, London, UK.

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, CLASP, CLASP II, the CLASP logo, Edwards PASCAL, Edwards PASCAL Precision, PASCAL and PASCAL Precision are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2025 Edwards Lifesciences Corporation. All rights reserved. PP--EU-9903 v1.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com



Edwards