



Sterilization & Disinfection

KynarPure-D can maintain integrity after being exposed to a variety of harsh sterilization methods, for results see table below.

Autoclave	124 °C, 30 minutes, 30 Cycles
Autoclave	134 °C, 30 minutes, 20 Cycles
Steam in Place	124 °C, 30 minutes, 30 Cycles
Steam in Place	134 °C, 30 minutes, 20 Cycles
Hot Water	85 °C, 60 minutes, 10 Cycles

The common sterilization and disinfection methods above have been tested under controlled laboratory conditions.

Note: This Cartridge has a PVDF membrane pre-filter laid over the PVDF membrane final filter and both layers are pleated together.

