


FLEXcon® OMNI-WAVE™ is a hydrogel-free biosensing technology.




- OMNI-WAVE™ passes ISO 10993-10 for sensitization and irritation.
- AAMI compliant.
- It has been awarded more than 90 patents worldwide and is ready to change the electrode paradigm.

FLEXcon Healthcare

New Possibilities in Wearables, Diagnostics, and Wellness.



FEATURES



OMNI-WAVE™



HYDROGEL

Performance and Stability	Not affected by humidity	Affected by humidity
Functional Life	3-day wear	2-3 day wear, skin and activity dependent
Out-of-pack Life	6-month out-of-pack life	30 to 60-day out-of-pack life
Allergic Reaction	No hydrogel related adverse reaction to skin	No standard options for those patients who are allergic to hydrogel
Ease of Processing	Easy to run - dry 2-mil (50 µm) adhesive	Roll to roll 30-40 mil (750 µm - 1,000 µm) gel processing can be challenging
Activity	Shower/bathe, jog	Failure in water
Storage Conditions	0° to 85°F -17° to 30°C	50° to 90°F 10° to 32°C
Packaging	No expensive barrier packaging required	Moisture barrier packaging required (foil + film + paper) packaging required
Cost Savings	\$30/1000 electrodes	