


Fogarty® Clot Management Catheters

Reference guide




Fogarty Corkscrew catheter

- For extracting thrombus that is too resistant for an elastomeric balloon
- Used to dislodge adherent material in native vessels and synthetic grafts
- Corkscrew design increases the surface area for entrapping fibrous material
- Radiopaque and flexible

Fogarty Corkscrew catheter

Catheter	End cap label	Model	Size	Graft / artery match
Corkscrew catheter	ACC 4/6	140806	4F x 6 mm x 80 cm	Tapered end of graft; 4 mm - 5 mm graft; peripheral arteries
Corkscrew catheter	ACC 5/8	140808	5F x 8 mm x 80 cm	6 mm - 7 mm graft; peripheral arteries
Corkscrew catheter	ACC 6/10	1408010	6F x 10 mm x 80 cm	6 mm - 10 mm graft; peripheral arteries



Fogarty graft thrombectomy catheter

- Designed for removing organised and adherent thrombus from synthetic grafts
- Provides a larger diameter range than the Fogarty Corkscrew catheter
- The pulling force can be varied by adjusting the diameter

Fogarty graft thrombectomy catheter

Catheter	End cap label	Model	Size	Graft match
Graft thrombectomy catheter	GT 5	160245F	5F x 16 mm x 50 cm	Polyester woven grafts and/ or ePTFE 6 mm or greater
Graft thrombectomy catheter	GT 6	160246F	6F x 18 mm x 50 cm	Polyester woven grafts and/ or ePTFE 8 mm or greater



Edwards

Fogarty clot management catheters

Model	Usable length (cm)	Catheter Size (F)	Membrane/spiral diameter max. extended (mm)	Membrane/spiral diameter max. retracted (mm)
Fogarty Corkscrew catheters				
140806	80	4	2.7	6
140808	80	5	3.1	8
1408010	80	6	3.3	10
Fogarty graft thrombectomy catheters				
160245F	50	5	5.0	16
160246F	50	6	6.0	18

Model	Usable length (cm)	Catheter size (F)	Inflated balloon diameter (mm)	Max. size of deflated balloon (F)	Usable length (cm)	Max. liquid capacity (ml)	Max. recommended pull force (lbs.)	Compatible introducer (F)	Compatible guidewire (in.)
Fogarty arterial embolectomy catheters (tube pack)									
120602F	60	2	4	3.4	0.2	–	0.5	3.5	–
120403F	40	3	5	4.3	0.6	0.20	0.7	4	–
120803F	80	3	5	4.3	0.6	0.20	0.7	4	–
120404F	40	4	9	5.0	1.7	0.75	1.5	5	–
120804F	80	4	9	5.0	1.7	0.75	1.5	5	–
120805F	80	5	11	6.0	3.0	1.50	2.0	6	–
120806F	80	6	13	7.0	4.5	2.00	2.5	7	–
120807F	80	7	14	8.0	5.0	2.50	3.5	8	–
Fogarty Fortis arterial embolectomy catheters									
120404F	40	4	9	5.0	1.7	0.75	–	5	–
120804F	80	4	9	5.0	1.7	0.75	–	5	–
Over-the-wire Fogarty thru-lumen embolectomy catheters									
12TLW403F	40	3	5	4.5	0.4	0.15	0.7	4.5	0.018
12TLW803F	80	3	5	4.5	0.4	0.15	0.7	4.5	0.018
12TLW1203F	120	3	5	4.5	0.4	0.15	0.7	4.5	0.018
12TLW404F	40	4	9	5.0	1.2	0.50	1.5	5	0.025
12TLW804F	80	4	9	5.0	1.2	0.50	1.5	5	0.025
12TLW405F35	40	5.5	11	6.0	1.7	0.90	2.0	6	0.035
12TLW805F35	80	5.5	11	6.0	1.7	0.90	2.0	6	0.035
12TLW406F	40	6	13	7.0	2.3	1.40	2.5	7	0.035
12TLW806F	80	6	13	7.0	2.3	1.40	2.5	7	0.035
12TLW807F	80	7	14	8.0	2.8	1.60	3.5	8	0.038
Fogarty venous thrombectomy catheters									
320806F	80	6	12	7.7	3.0	1.50	–	–	–
320808F	80	8	13	8.8	4.5	2.25	–	–	–
32080810F	80	8	19	10.0	8.0	4.00	–	–	–

CAUTION: For professional use only. See Instructions for Use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards, Edwards Lifesciences, the stylised E logo, Fogarty, Fogarty Corkscrew, Fogarty Fortis, and 1-800-FOGARTY are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

©2022 Edwards Lifesciences Corporation. All rights reserved. ANZ-2022-065

Edwards Lifesciences Pty Ltd • 2/40 Talavera Road, North Ryde NSW 2113 Australia phone: 1800 222 601

Edwards Lifesciences (New Zealand) Ltd • PO BOX 28658 Remuera New Zealand • New Zealand Phone: 0800 222 601

