

OptiSite


Aortic Cannulae

Designed to go where you need it


Every patient and surgical situation is different, which is why you need an arterial cannula that provides options. OptiSite arterial cannulae incorporate an easy-to-use thin-wall technology to provide you with a variety of cannulation site options, designed to accommodate multiple surgical techniques and clinician preferences. Indicated for cannulation in various arterial sites such as the aorta, femoral, axillary and subclavian.

- Smooth rounded introducer tip allows for atraumatic insertion
- Locking introducer reduces retraction during insertion
- Vent cap reduces back-bleed during aortic cannulation
- Depth markings aid in placement
- Fr. sizes include 16, 18, 20 and 22

Pressure Drop vs Flow



Smooth Rounded Tip of the OptiSite arterial cannula



Vent Cap reduces back-bleed



OptiSite cannulae models

Model Number	Size	Effective Length	Connector
OPTI16	16 Fr. (5.3 mm)	6.0" (15 cm)	3/8" Vented
OPTI18	18 Fr. (6.0 mm)	6.0" (15 cm)	3/8" Vented
OPTI20	20 Fr. (6.7 mm)	6.0" (15 cm)	3/8" Vented
OPTI2	22 Fr. (7.3 mm)	6.0" (15 cm)	3/8" Vented



Insertion for additional arterial access sites such as the aorta, axillary, and subclavian

Insertion Sites	OptiSite
Axillary	✓
Aorta	✓
Subclavian	✓
Femoral	✓

CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards, Edwards Lifesciences, the stylized E logo, and OptiSite are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

© 2021 Edwards Lifesciences Corporation. All rights reserved. PP--US-1555 v3.0

Edwards Lifesciences • One Edwards Way, Irvine CA 92614 USA • edwards.com

