


Edwards PASCAL Precision Transcatheter Valve Repair System









PASCAL and PASCAL Ace Implant System

Technical Specifications

Disposable - single use (sterile)

Implant System		
		
Components	Reference number	
1 PASCAL implant	20000IS	
1 PASCAL Ace implant	20000ISM	
2 Steerable catheter	3 Implant catheter	

1 PASCAL Implant / PASCAL Ace Implant

	Closed	Leaflet - Capture - Ready		Elongated
PASCAL Ace				
PASCAL				
	Medial-Lateral	Medial-Lateral	Anterior-Posterior	Anterior-Posterior

Note: Implant dimensions are in millimeters and are approximate for reference.


Disposable - single use (sterile)

Accessory Box

1 Guide Sheath Set (sterile)

	Components	Reference number
	1. Guide sheath 2. Introducer	20000GS

2 Stabilizer (sterile)





	Reference number
	20000STP

System Specifications

Component	Dimension
Implant Catheter	
Implant catheter working length ^A	118 cm
Steerable Catheter	
Steerable catheter working length ^B	109 cm
Steerable catheter outer diameter (OD)	5.7 mm
Guide Sheath and Introducer	
Guide sheath working length	79 cm
Guide sheath proximal shaft outer diameter (OD)	7.3 mm / 22 Fr
Guide sheath proximal shaft inner diameter (ID)	5.8 mm
Introducer working length	104 cm
Introducer inner diameter (ID)	1.0 mm (0.035" Guidewire compatible)

Reusable Accessories

Accessory Box

1 Reusable Cradle (non-sterile)		2 Reusable Plate (non-sterile)	
			
Reference number		Reference number	
10000CR		10000PT	
3 Reusable Platform (non-sterile)			
Legs Not Extended		Legs Fully Extended	
		 Maximum height is 28.3 cm	
Legs Retracted for Storage		Reference number	
		10000UP	

Note: Since these are reusable, they are not required to be reordered for every procedure.

PASCAL Precision System - Single Use

	Implant type	Customer order number	Devices	Reference number
Edwards PASCAL Precision Transcatheter Valve Repair System with PASCAL Implant	PAS2-M10* PAS2-T10*	20000IS	Implant Box	
			Implant system	20000IS
	PAS2-M10* PAS2-T10*	21000GSKIT	Accessory Box	
			1 x Guide sheath	20000GS
			1 x Stabilizer	20000STP
Edwards PASCAL Precision Transcatheter Valve Repair System with PASCAL Ace Implant	PAS2-M5* PAS2-T5*	20000ISM	Implant Box	
			Implant system	20000ISM
	PAS2-M5* PAS2-T5*	21000GSKIT	Accessory Box	
			1 x Guide sheath	20000GS
			1 x Stabilizer	20000STP

*M10 or M5 for mitral procedures, T10 or T5 for tricuspid procedures.

Reusable Accessories

	Note	Customer order number	Devices	Reference number
Reusable Accessories	Order this code to supply your account with the Reusable Accessories Starter Kit, including 2x of each component.	10000RAKIT	Reusable Accessory Starter Kit	
			2 x Reusable Cradle	10000CR
			2 x Reusable Plate	10000PT
			2 x Reusable Platform	10000UP
	If replacement individual Reusable Accessories are required, order the individual component code to restock the item.			

If you have any questions, please contact your local Edwards sales representative.

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, Edwards PASCAL, Edwards PASCAL Precision, PASCAL, PASCAL Ace and PASCAL Precision are trademarks or service marks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

© 2023 Edwards Lifesciences Corporation. All rights reserved. PP--EU-6117 v1.0

Edwards Lifesciences • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com



Edwards