


SAPIEN M3


Transcatheter Mitral Valve Replacement System

Technical Specifications


SAPIEN M3 System



1 SAPIEN M3 Dock



2 SAPIEN M3 Valve






Reference
Number

N/A
(See SAPIEN M3 Dock
Steerable Catheter)

9880TFX29MCE

Delivery System(s) - Disposable - Single Use (Sterile)

1 Edwards 23F Guide Sheath	Components	Reference Number
	1a Edwards 23F Guide Sheath	9880GSCE
	1b Introducer	
2 SAPIEN M3 Dock Steerable Catheter (includes SAPIEN M3 Dock)	Components	Reference Number
	SAPIEN M3 Dock Steerable Catheter	9880DDSCE
	SAPIEN M3 Dock	
	Preparation Accessories	
3 Edwards Commander M Delivery System	Components	Reference Number
	Edwards Commander M Delivery System	9880CM29CE
	Loader	
	Crimp Stopper	


	Working Length	ID	OD
1a Edwards 23F Guide Sheath	77 cm	23F	29F
1b Introducer	103 cm	N/A	24F
2 SAPIEN M3 Dock Steerable Catheter	113 cm	N/A	18F
3 Commander M Delivery System	112 cm	N/A	16F



Disposable Accessories - Single Use (Sterile)

SAPIEN Stabilizer Rail System		SAPIEN M3 Crimper	Edwards Locking Syringe (Inflation Device)
			
Reference Number	9880SR5CE	9880CRCE	96406

Reusable Accessories - Non-sterile

Note: Since these are reusable, they are not required to be reordered for every procedure.

Edwards Reusable Platform (Adjustable)	
Legs Not Extended	Legs Fully Extended
	
Reference Number	10000UP

Edwards Reusable Plate	Edwards Reusable Cradle
	
Reference Number	10000CR

SAPIEN M3 System Ordering Information - Single Use

Customer Order Number	Box Description	Components	Reference Number
SM329EU1	SAPIEN M3 Transcatheter Mitral Valve Replacement System	Edwards 23F Guide Sheath	9880GSCE  (01)00690103218138
		SAPIEN M3 Dock Steerable Catheter (includes SAPIEN M3 Dock)	9880DDSCCE  (01)00690103218114
		SAPIEN M3 Valve	9880TFX29MCE  (01)00690103218107
		Edwards Commander M Delivery System	9880CM29CE  (01)00690103218121
		SAPIEN Stabilizer Rail System	9880SRSCCE  (01)00690103218152
		SAPIEN M3 Crimper	9880CRCE  (01)00690103218145
		Edwards Locking Syringe (Inflation Device)	96406  (01)00690103210019

Reusable Accessories Ordering Information - Non-sterile

Note: Since these are reusable, they are not required to be reordered for every procedure. If replacing an individual reusable accessory is required, order the individual reference number to restock the item.

Customer Order Number	Box Description	Components	Reference Number
10000RAKIT*	Reusable Accessory Starter Kit*	Edwards Reusable Cradle	10000CR  (01)00690103218015
		Edwards Reusable Plate	10000PT  (01)00690103217278
		Edwards Reusable Platform	10000UP  (01)00690103217292

*The Reusable Accessory Starter Kit includes two of each component.

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, Commander, Commander M, Edwards Commander, Edwards Commander M, Edwards SAPIEN, Edwards SAPIEN M3, SAPIEN, and SAPIEN M3 are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2025 Edwards Lifesciences Corporation. All rights reserved. PP--EU-9235 v2.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com



Edwards