


# **COVID-19 Brady Workstation Software label layouts**

Attached you will find a list of free COVID-19 related Brady Workstation Software label layouts and their corresponding compatible Brady benchtop printer and material width.

Download **COVID-19 Brady Workstation Software label layouts** here: <u>BradyID.com/COVID19</u>

For additional support please contact Brady Technical Support: 1-800-643-8766



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 4 in. Continuous COVID-19 Notice Guidelines	4 in.



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP(35)37 4in x 6in COVID-19 Clean High Touch Areas	4 in.



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP(35)37 4in x 6in COVID-19 Clean Surfaces Notice	4 in.
BradyPrinter S300 and S3100	S3100(S3000) 2 in Continuous COVID-19 Clean Surfaces Notice	2 in.



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 4in x 6in COVID-19 Slow the Spread	4 in.
BBP®85	BBP85 10in Continuous COVID-19 Slow the Spread	10 in.
BradyJet J5000	J5000 8in x 11in COVID-19 Slow the Spread	8 in.
BradyPrinter S3100 and S3000	S3100(S3000) 2in Continuous COVID-19 Slow the Spread	2 in.



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP(35)37 4in x 6in COVID-19 Keep Hands Clean	4 in.
BradyJet J5000	J5000 8in x 11in COVID-19 Slow the Spread	8 in.
BradyPrinter S3100 and S3000	S3100(S3000) 2in Continuous COVID-19 Slow the Spread	2 in.




Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 4in Continuous COVID-19 Attention Guidelines	4 in.

## Please Use Hand Sanitizer






Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP(35)37 4in x 6in COVID-19 Hand Sanitizer Notice	4 in.
BradyPrinter S3100 and S3000	S3100(S3000) 2in COVID-19 Hand Sanitizer Notice	2 in.



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 4in x 6in COVID-19 Social Distance Notice	4 in.
BBP®85	BBP85 10in Continuous COVID-19 Social Distance Notice	10 in.
BradyJet J5000	J5000 8in x 10in COVID-19 Social Distance Notice	8 in.



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 4in x 6in COVID-19 Symptom Checker	4 in.
BradyJet J5000	J5000 8in x 11in COVID-19 Symptom Checker	8 in.
BradyPrinter S3100 and S3000	S3100(S3000) 2in Continuous COVID-19 Symptom Checker	2 in.



Printer	Brady Workstation file name	Material width
BradyJet J5000	J5000 4in COVID-19 Wash Hands Notice	8 in.



Printer(s)	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 4in Continuous COVID-19 Safety First Social Distance	4 in.




Printer	Brady Workstation file name	Material width
BradyJet J5000	J5000 8in COVID-19 Awareness Sign	8 in.



Printer	Brady Workstation file name	Material width
BradyJet J5000	J5000 8in x 11in COVID-19 Clean High Touch Areas	8 in.

#### Floor tape:



Printer	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 2.25in COVID-19 Floor Tape 6ft Apart	2.25 in.




Printer	Brady Workstation file name	Material width
BBP®35 BBP®37	BBP35(37) 2.25in Floor Tape Please Wait	2.25 in.

### How to download COVID-19 Brady Workstation Software label layouts from BradyID.com


- 1. Go to <a href="mailto:BradyID.com/COVID19">BradyID.com/COVID19</a>.
- 2. Download the folder to your Desktop.
- 3. Right click on the folder COVID-19 Workstation Label Layouts, choose Extract All.
- Double click on extracted COVID-19 Workstation Label Layouts folder.
- 5. Select the folder with your printer.
- 6. Double click on the Brady Workstation label file you need.

#### How to change the label part or printer in Brady Workstation

1. In the Document Properties Tab, click Change Part or Printer.



2. Click the dropdown menu and select your printer.



#### How to change the label part or printer in Brady Workstation

- 3. Select your label part after selecting your printer.
- 4. Click save to confirm your changes and go back to the design page to create / modify your label.

