



Medical Director printer setting instructions for I-MED Radiology

1. Click **Tools > Print options > Imaging**.
2. Select Provider **I-MED Radiology** from the drop down list.
3. Tick **Use as default**.
4. Suggested Font/size = **Arial 10**.
5. Follow the layout below and click **Save** to update print options.

	From left	From top		From left	From top
Patient's Name	70	48	Doctor's Name	80	240
Patient's Address	70	53	Dr's address	80	245
Patient's DOB	170	48	Dr's Provider no.	165	240
Medicare No.	170	61	Dr's Phone	80	258
Chart Number	170	67	Date Printed	160	265
Phone No.	70	61	Copy to	165	253
	From left	From top	Width	Lines	
Request	70	110	120	3	
Reason	60	140	140	14	

Reset defaults

All measurements are in millimeters from the top left corner of the page

Display default reason list?

Save

Cancel

Help

If you have any troubles, please reach out to your I-MED Radiology Account Manager or contact our Referrer Success Team. Visit <https://i-med.com.au/support>