



# Medical Director printer setting instructions for I-MED Radiology

1. Click **Tools > Print options > Imaging**.
2. Select Provider **I-MED Radiology** from the drop down list.
3. Tick **Use as default**.
4. Suggested Font/size = **Arial 10**.
5. Follow the layout below and click **Save** to update print options.

	From left	From top		From left	From top
Patient's Name	35	25	Doctor's Name	35	230
Patient's Address	35	45	Dr's address	35	240
Patient's DOB	35	40	Dr's Provider no.	85	230
Patient's Sex	85	25			
Medicare No.	105	25	Dr's Phone	85	240
Chart Number	0	0	Date Printed	140	260
Phone No.	35	60	Copy to	35	255
	From left	From top	Width	Lines	
Request	35	115	150	50	
Reason	35	165	150	60	

Reset defaults

All measurements are in millimeters from the top left corner of the page

Display default reason list?

Save

Cancel

## Help

If you have any troubles, please reach out to your I-MED Radiology Account Manager or contact our Referrer Success Team. Visit <https://i-med.com.au/support>