



I-MED Online install instructions for MAC

Before commencing install:

- Check that the Apple operating system is Version 10.6.0 or greater
- Upgrade software if version is less

Choosing a default web browser – Safari is recommended:

- Select Apple menu (top left of screen)
- In that menu select 'System Preferences'
- Go into 'General'
- Choose your desired browser from the 'Default Web Browser' drop down menu

How to open shortcut in Safari when the default browser is something else:

- Ctrl-click a shortcut (or select more than one) & ctrl-click, choose 'Get Info'
- In the 'Open With' section, change it to Safari
- DO NOT click the 'Change All' button just close the info window(s). That should force just those particular shortcuts to use Safari

Click on the below link or type the following into your web browser: https://i-med.com.au/doctors/access-patient-images

On this page select Regional NSW, then the installers button.

- 1. Click on the MAC installer to download the Inteleviewer.dmg file.
- 2. Locate the file in your downloads and double click on the file to commence the install.
- 3. Open a new finder window (right click Finder), and select the applications folder.
- 4. Click and drag the Inteleviewer icon from step 2, into the applications folder.
- 5. Close all open files, then open safari and type: <u>https://i-medonline.ril.com.au/webcomrad/</u> to locate the I-MED online website.
- 6. Click on bookmarks and save as a bookmark, then rename (webCOMRAD) I-MED Online.
- 7. Drag the I-MED Online bookmark to the desktop.
- 8. Click to open I-MED Online (using Safari) and Log in using your secure User ID and Password.

For additional support please contact your Customer Relationship Manager at I-MED Radiology.



MAC Installer (ZIP 286 MB)





I-MED Online install instructions for MAC (Cont.)

I-MED Online has been updated to resolve the ongoing issues referring practitioners have had whenever Java updates were released by Apple, as this was preventing InteleViewer from launching.

There are a few steps that need to be followed to install a launcher application on a referring practitioner's Mac, and then InteleViewer should launch, regardless of Java updates.

Step 1: Log in to I-MED Online and view a report

Step 2: Click on the 'Download Inteleviewer Launcher for Mac' link just underneath the 'Online Images' button



Step 3: Once the download has finished, double-click on the 'ComradInteleviewer.dmg' file that has been downloaded

Step 4: In the window that opens, drag the 'Inteleviewer Launcher' file onto 'Applications'



Step 5: Click on the 'Launch InteleViewer' button in I-MED Online. If the InteleViewer loads, you're finished! If instead you see this error:



Keep following these instructions (Continue to next page).





I-MED Online install instructions for MAC (Cont.)

Step 6: Open a new Finder window

Finder	File	Edit	View	Go	Window	Help	
	New Finder Window New Folder						ℋN
							0 X N
	New Folder with Selection					~ 36 N	
	New Smart Folder						X 96 N

Step 7: Click on Applications in the Finder window, and find the InteleviewerLauncher file

	Applica	tions		
4 >	# E III # & # # #	Q	X	je i
FAVORITES	Name	Date Modified	Size	Ki
All My Files	@ GarageBand	12/07/2011 4:19 AM	381.7 MB	AL
C AirDonn	觉 iCal	01/10/2012 3:34 PM	56.8 MB	Ap
T Minorop	O iChat	01/10/2012 3:34 PM	51 MB	At
Applications	Image Capture	01/10/2012 3:34 PM	5.5 MB	Ap
Desktop	yk iMovie	12/07/2011 4:31 AM	695.3 MB	Ap
P Documents	11 InteleViewer	10/09/2011 6:15 AM	25 MB	Ap
O Downloade	C InteleviewerLauncher	12/08/2013 9:09 AM	103 KB	Ap
Downloads	🚮 iPhoto	12/07/2011 4:32 AM	1.09 GB	Ap
Movies	A	1010010010 0.17 111	200 642	

Step 8: Right click on InteleviewerLauncher, and click Open. You'll see the following warning:



Step 9: Click Open button. Another warning will appear:



Step 10: Click the OK button

Step 11: Quit Safari

Now click on the 'Launch Inteleviewer' button in I-MED Online – Inteleviewer should open successfully!

For additional support please contact your Customer Relationship Manager at I-MED Radiology.