



Medical Director printer setting instructions for I-MED Radiology

1. Click **Tools > Print options > Imaging**.
2. Select Provider **I-MED Radiology** from the drop down list.
3. Tick **Use as default**.
4. Suggested Font/size = **Arial 10**.
5. Follow the layout below and click **Save** to update print options.

	From left	From top		From left	From top
Patient's Name	30	37	Doctor's Name	30	105
Patient's Address	30	49	Dr's address	30	110
Patient's DOB	170	37	Dr's Provider no.	30	122
Medicare No.	170	53	Dr's Phone	60	122
Chart Number	170	65	Dr's Fax	60	122
Phone No.	170	43	Date Printed	60	27
My Health Record	118	195	Copy to	145	120
eOrder barcode	0	0			
	From left	From top	Width	Lines	
Request	30	82	75	6	
Reason	115	82	75	6	

Reset defaults

All measurements are in millimeters from the top left corner of the page

Display default reason list?

Save Cancel

Help

If you have any troubles, please reach out to your I-MED Radiology Account Manager or contact our Referrer Success Team. Visit <https://i-med.com.au/support>